

Advanced Vibration Technology to Improve Patient Wellbeing

pro7HC™



As the leading brand of whole body vibration equipment, Power Plate has been scientifically proven to address a range of medical conditions, speed recovery time, and quickly return patients to their activities of daily living.

Our newest model, the pro7 Healthcare, is the first commercial vibration training machine designed specifically for the needs of medical facilities and treatment centers. The pro7HC comes complete with patient friendly, integrated LCD touch screen video programs focused on neurological (neuro) health, pain relief and common conditions. Equipped with an integrated heart rate monitor and embedded multidirectional cables, range of motion tracking, and real-time variable resistance (power output) display, the pro7HC makes assisting patients easier and more efficient.

Whether you are looking to address neurological based conditions, physical pain, or general health issues, the pro7HC is a valuable addition to any practice.



Power Plate® pro7HC™ Technical Specifications

71-HC7-3300

Frequency	25-50Hz
Amplitude	Low or High
Vibration Settings	0–8 (with 27 sublevels of intensity)
Time Selections	0–9 minutes (15 second increments)
DualSync™Twin Motor System	Maintains precise balance at any frequency and amplitude level, allowing perfect synchronization of vibration for maximum muscle response and efficiency
Total Dimensions (W x D x H)	38in x 46in x 61in / 96cm x 116cm x 155cm
Platform Dimensions (W x D x H)	38in x 37in / 96cm x 95cm x 30cm
Weight	416lb / 189 kg
Power Supply	100–240V, 50/60Hz, Nominal Power: 160-265W, Universal Voltage
Maximum Load	500lb / 227kg
proMOTION™ Dynamic Vibration Technology	Touch screen tension adjustment (three levels) with real-time resistance display; embedded with modular attachment option
Stabilization System	Advanced Stabilization and Compensation System
Heart Rate Monitor	ANT+ wireless technology
Operating System	Windows CE 6.0 Professional
Upper Display	10.1in LCD Touchscreen Display (800 x 480 dpi)
USB Port/Hub	USB 1.0 Port/ 1 external, 1 internal service port
Net Access	10/100M Ethernet
Audio	2 recessed speakers with headphone option and volume control
Certifications	MDD, CE and EMC; RoHS/WEEE, NRTL, C-NRTL, FDA listed as Class 1 device, 510k exempt, CB
Warranty	2 years parts 1 year labor

Performance Health Systems

401 Huehl Rd., Suite 2A
 Northbrook, IL 60062
 (877) 877-5283 | info@powerplate.com | powerplate.com